

HIGH-PRECISION ION MOBILITY CALIBRATION IMPROVES DIA-PASEF ANALYSIS



Tejas Gandhi, PhD
Bioinformatics Manager

tejas.gandhi@biognosys.com
www.biognosys.com

Tejas Gandhi¹, Stephanie-Kaspar Schoenefeld², Oliver M. Bernhardt¹, David Schlessinger¹, Sven Brehmer², Gary Krappa², Jan Muntel¹, Lynn Verbeke¹, Lukas Reiter¹


1) Biognosys AG, Wagistrasse 21, 8952 Schlieren (Zurich), Switzerland 2) Bruker Daltonik GmbH, Bremen, Germany

INTRODUCTION

Ion mobility (IM) based data independent acquisition has recently become available on timsTOF Pro. Due to its parallel accumulation-serial fragmentation (PASEF) strategy it achieves a high duty cycle and has shown promise in biological applications. We have previously shown that high-precision retention time calibration benefits DIA analysis (Bruderer, 2016). Upon applying this concept

to ion mobility for analyzing dia-PASEF data, we could identify a high number of proteins, competitive to state-of-the-art DIA methods.

Figure 1: High-precision IM in Spectronaut
Ion mobility calibration (left panel) ensures that optimal IM extraction width (red line, right panel) is automatically used for data analysis as shown here with a 400 ng HeLa sample.



CONCLUSIONS

- A novel workflow for high-precision ion mobility calibration in Spectronaut improves identification (Figure 1) and quantification (Figure 3)
- Deep protein coverage achieved even with low sample amounts (Figure 2)
- SN 14 improves on its dia-PASEF workflow from SN 13 (Figure 3)
- Spectronaut processes timsTOF Pro data significantly faster than acquisition time (Figure 2)

RESULTS




Figure 2: Triplicate HeLa lysate dia-PASEF runs with varying sample amounts (10ng, 50ng, 400ng) on a 100 min gradient

Top Panel: Average number of proteins identified per injected sample amount with 1% peptide and protein FDR. The library consisting of over 10,000 protein groups was created from 15 fractions using Pulsar search engine with 1% FDR (PSM, precursor and protein)

Bottom Panel: The analysis time per run (in minutes) for the library generation (PASEF) and subsequent, targeted analysis of the 400 ng dia-PASEF runs is significantly lower than the acquisition time (100 min gradient + 20 min overhead)




Figure 3: UPS2 proteins spiked in 50 ng mouse cerebellum background matrix in 5 different steps, 2 replicates each

Top Panel: Average number of proteins identified by different Spectronaut versions. The red line indicates the library size in proteins.

Bottom Panel: Regulation analysis was performed using an unpaired t-test of protein quantities calculated from both MS1 and MS2 (Huang, 2019) in SN 14. For SN 13, the t-test was based on precursor quantities calculated from MS2. Upon sorting all the regulated candidate pairs by p-value, we found 34 UPS2 candidate pairs in the top 100. This is comparable to what we have previously observed in this sample (Muntel, 2019).