

Untargeted UHPLC-HR-QTOF-MS metabolomics study to unravel metabolites controlling wood formation in aspen trees



ASMS 2020, WP 444

Magdalene Reinkensmeier¹, Ilara G.F. Budzinski², Thomas Moritz^{3,4}

¹Bruker Daltonics, Bremen, Germany
²Umeå Plant Science Centre, Department of Plant Physiology, Swedish University of Agricultural Sciences, Umeå, Sweden
³Swedish Metabolomics Centre, Swedish University of Agricultural Sciences, Umeå, Sweden
⁴Novo Nordisk Foundation Centre for Basic Metabolic Research, University of Copenhagen, Copenhagen, Denmark


Introduction

Plants as sessile organisms have to cope with the environmental conditions in their inhabited soil, e.g. the existing amount of water and nutrients. The available resources must be allocated to growth and the handling of biotic challenges, such as the defense against predators. Secondary growth in trees is of great economic importance as it results in the production of wood. The developmental process of wood takes place in the vascular cambium. Cambium activity results in the regular production of xylem and phloem, in opposite directions, where xylem having a role in water transport and tree stability, phloem has a role in sugar transport. Wood formation is a dynamic and continuous process, strongly affected by environmental factors such as water, nutrient, and light availability. We investigated the impact of variations in nitrogen fertilization on the growth and the metabolite profiles of xylem and phloem of aspen trees using untargeted LC-HR-QTOF-MS metabolomics. Within the untargeted approach we have put a focus on phenolic glucosides, the salicylates, which are major defense secondary metabolites of aspen and other *Salicaceae* species.


Methods

Cuttings of hybrid aspen (*Populus tremula* x *P. tremuloides*) were cultivated in a greenhouse. After two months the trees were treated with low or high nitrogen fertilizer three times per week for two months. At the harvest stem segments were sampled, flash-frozen in liquid nitrogen, divided into xylem and phloem samples with scraping, and thereafter stored at -80°C. Samples (10 mg each) were extracted according to Gullberg et al. 2004 [1]. Samples were analyzed with the Elute UHPLC (ACQUITY UPLC HSS T3, 2.1 x 50 mm, 1.8 µm, C18 column and 2.1 mm x 5 mm, 1.8 µm, VanGuard Pre-column (Waters, MA, USA) coupled to the impact II HR-QTOF-MS (Bruker Daltonics, MA, USA). Default metabolomics acquisition methods were used to collect full scan and fragment data in electrospray positive and negative ion mode. MetaboScape[®] 5.0 was used for comprehensive feature finding with the T-ReX[®] 3D algorithm, feature annotation, statistical data evaluation, and pathway mapping. Analyte lists created from in-house databases were used to annotate target compounds (including e.g. phenolic glucosides, oligolignols, amino acids). Pathway maps were prepared in PathVisio 3.3.0. Additional measurements with the timsTOF Pro (Bruker Daltonics, MA, USA) were done to collect ion mobility data.


Workflow



Results





→ Considerable differences in the metabolite profiles




PCA scores (left) and loadings plot (right) for the merged metabolite profiling data ($N=3$ per group). Annotated salicylates and related compounds are named. Salicylates shown in the pathway maps are in bold.

→ A closer look on salicylates




→ Pathway mapping illustrate cell type specific metabolic response



→ Enhanced annotation certainty through inclusion of CCS values

The addition of CCS value into the annotation workflow in MetaboScape[®] led to an elevated annotation certainty for numerous features. Reference CCS values were extracted from the Unified CCS Compendium [6].

Molecular For...	Name	AQ	CCS (Å ²)	ΔCCS [%]
C ₅ H ₅ N ₅	Adenine	■	124.1	1.2
C ₉ H ₉ N ₂ O ₂	L-Phenylalanine	■	138.9	1.4
C ₁₀ H ₁₃ N ₅ O ₄	Adenosine	■	154.9	1.8
C ₁₅ H ₁₄ O ₆	Catechin	■	168.0	2.4



References

- Gullberg J. et al., Anal. Biochem. 2004, 331:283-295
- Wolf S. et al., BMC Bioinformatics 2010, 2010:1148
- Ruttikies C. et al., J. Cheminformatics 2016, 8:3
- Abreu I.N. et al., J. Chem. Ecol. 37(8):857-870
- Snyder D.T. et al., Anal. Methods 2015, 7:870-876
- Picache, J. A. et al., Chem. Sci. 2019, 10:983-993

Conclusions

- Different amounts of N fertilizer led to distinct changes in metabolite profiles for both phloem and xylem tissues.
- Salicylates were strongly affected especially in xylem.
- Other metabolite groups such as oligolignols also showed specific patterns in the different sample cell types and fertilizer groups (data not shown here).
- This pilot study shows that the experimental set-up, in combination with untargeted LC-MS metabolomics analysis, has potential to identify metabolite markers controlling important aspects of plant growth and development.
- Further experiments with TIMS technology will help to get a deeper insight into the complex xylem and phloem samples and lead to more comprehensive feature annotation by incorporating a further dimension of resolution.

UHR-QqTOF