

From hair to body liquids – TXRF analysis of biological and medical samples



Bruker Nano Analytics, Berlin
Webinar, December 18th, 2019



Welcome




Speakers

Dr. Hagen Stosnach
Applications Scientist TXRF
Berlin, Germany



Prof. Dr. Lutz Schomburg Professor for
Experimental Endocrinology and Deputy Director
Charité-University medicine Berlin, Institute for
experimental endocrinology



Part I: Analysis of human hair

- Motivation
- TXRF spectroscopy

Part II:

- Bone as matrix
- Se-Supplementation Trial in Humans
- Diagnostic Trial of serum Se and Cu in TSCI
- Analytics of recombinant or purified Protein
- Nested case-control Study => predictive Biomarker

Part III: Summary/Q & A



Part I: Analysis of human hair

Analysis of human hair

Motivation



Significance of micronutrients in human hair

The distribution of micronutrients like Mg, K, Ca, Mn, Cu, Zn, Se, Mo and I gives information on diseases, metabolic disorders, environmental exposures, and nutritional status.

Examples:

- Breast cancer patients show low hair Ca, Mg and Zn, but high As, Na and K concentrations compared with the normal controls. The Fe level is shown to be significantly low and associated with Ca and Mn levels [1]
- A cross-sectional analysis on 343 subjects showed that those with metabolic syndrome had significantly lower contents of Ca, Mg and Cu, whereas the amounts of Na, K and Hg are higher [2]
- Hair Zn, Se and Cu deficiencies are noted in chronic gastrointestinal diseases [3]

Analysis of human hair

Motivation



Control of toxic metals in human hair

- Hair analysis has been considered as one of the most important biomarkers according to the Environmental Protection Agency (EPA) [4].
- As an excretory system, human hair can accumulate and incorporate toxic metals into its structure during its growth process. Therefore, concentrations of heavy metals in hair can reflect the mean level in the human body, recording the population's exposure to heavy metals [5]

[1] Joo NS, Kim SM, Jung YS, Kim KM. Hair iron and other minerals' level in breast cancer patients. *Biol Trace Elem Res.* 2009;129:28–35.

[2] Park SB, Choi SW, Nam AY. Hair tissue mineral analysis and metabolic syndrome. *Biol Trace Elem Res.* 2009;130:218–28.

[3] Bhat YJ, Manzoor S, Khan AR, Qayoom S. Trace element levels in alopecia areata. *Indian J Dermatol Venereol Leprol.* 2009;75:29–31

[4] Rashed M.N., Hossam F. Heavy metals in fingernails and scalp hair of children, adults and workers from environmentally exposed areas at Aswan. Egypt. *Environ. Bioindic.* 2007;2:131–145

[5] Sera K., Futatsugawa S., Murao S. Quantitative analysis of untreated hair samples for monitoring human exposure to heavy metals. *Nucl. Instrum. Methods B.* 2002;189:174–179.

Analysis of human hair

Motivation



Benefits of human hair analysis

- Hair is easier to collect, transport and store than other medical samples like blood or serum
- As elements are fixed in the hair, element control over longer time period is possible


Analyzed samples

- Certified reference sample NIES 13 (National Institute for Environmental Studies)
- 5 human hair samples (cleaned with acetone and ultrapure water)



Analysis of human hair

Samples and sample preparation



Analysis of human hair Measurements



S2 PICOFOX

- Mo tube, 50 kV/1000 μ A
- 60 mm² XFlash SDD

Measurement conditions

- Mo-K excitation, 1000 s



Analysis of human hair Measurements



Measurement program

1 H Hydrogen																	2 He Helium						
3 Li Lithium	4 Be Beryllium	<div style="display: flex; align-items: center; gap: 20px;"> <div style="background-color: #0070C0; width: 20px; height: 20px; display: inline-block;"></div> Analysed using K-lines <div style="background-color: #00B050; width: 20px; height: 20px; display: inline-block;"></div> Analysed using L-lines </div>										5 B Boron	6 C Carbon	7 N Nitrogen	8 O Oxygen	9 F Fluorine	10 Ne Neon						
11 Na Sodium	12 Mg Magnesium																	13 Al Aluminium	14 Si Silicon	15 P Phosphorus	16 S Sulphur	17 Cl Chlorine	18 Ar Argon
19 K Potassium	20 Ca Calcium	21 Sc Scandium	22 Ti Titanium	23 V Vanadium	24 Cr Chromium	25 Mn Manganese	26 Fe Iron	27 Co Cobalt	28 Ni Nickel	29 Cu Copper	30 Zn Zinc	31 Ga Gallium	32 Ge Germanium	33 As Arsenic	34 Se Selenium	35 Br Bromine	36 Kr Krypton						
37 Rb Rubidium	38 Sr Strontium	39 Y Yttrium	40 Zr Zirconium	41 Nb Niobium	42 Mo Molybdenum	43 Tc Technetium	44 Ru Ruthenium	45 Rh Rhodium	46 Pd Palladium	47 Ag Silver	48 Cd Cadmium	49 In Indium	50 Sb Antimony	51 Sn Tin	52 Te Tellurium	53 I Iodine	54 Xe Xenon						
55 Cs Cesium	56 Ba Barium	57 La Lanthanum	72 Hf Hafnium	73 Ta Tantalum	74 W Tungsten	75 Re Rhenium	76 Os Osmium	77 Ir Iridium	78 Pt Platinum	79 Au Gold	80 Hg Mercury	81 Tl Thallium	82 Pb Lead	83 Bi Bismuth	84 Po Polonium	85 At Astatine	86 Rn Radon						
87 Fr Francium	88 Ra Radium	89 Ac Actinium																					
			L Lanthanides	58 Ce Cerium	59 Pr Praseodymium	60 Nd Neodymium	61 Pm Promethium	62 Sm Samarium	63 Eu Europium	64 Gd Gadolinium	65 Tb Terbium	66 Dy Dysprosium	67 Ho Holmium	68 Er Erbium	69 Th Thulium	70 Yb Ytterbium	71 Lu Lutetium						
			Ac Actinides	90 Th Thorium	91 Pa Protactinium	92 U Uranium	93 Np Neptunium	94 Pu Plutonium	95 Am Americium	96 Cm Curium	97 Bk Berkelium	98 Cf Californium	99 Es Einsteinium	100 Fm Fermium	101 Md Mendelevium	102 No Nobelium	103 Lr Lawrencium						

Analysis of human hair Measurements



S4 T-STAR

- Mo tube, 50 kV/1000 μ A
- W-tube, 50 kV/1000 μ A
- 60 mm² XFlash SDD

Measurement conditions

- Mo-K excitation, 1000 s
- W-Brems excitation, 1000 s



Analysis of human hair Measurements




Excitation modes used for quantification

1																		2																	
H Hydrogen																		He Helium																	
3																		5		6		7		8		9		10							
Li Lithium																		B Boron		C Carbon		N Nitrogen		O Oxygen		F Fluorine		Ne Neon							
11																		13		14		15		16		17		18							
Na Sodium																		Al Aluminium		Si Silicon		P Phosphorus		S Sulphur		Cl Chlorine		Ar Argon							
19		20		21		22		23		24		25		26		27		28		29		30		31		32		33		34		35		36	
K Potassium		Ca Calcium		Sc Scandium		Ti Titanium		V Vanadium		Cr Chromium		Mn Manganese		Fe Iron		Co Cobalt		Ni Nickel		Cu Copper		Zn Zinc		Ga Gallium		Ge Germanium		As Arsenic		Se Selenium		Br Bromine		Kr Krypton	
37		38		39		40		41		42		43		44		45		46		47		48		49		50		51		52		53		54	
Rb Rubidium		Sr Strontium		Y Yttrium		Zr Zirconium		Nb Niobium		Mo Molybdenum		Tc Technetium		Ru Ruthenium		Rh Rhodium		Pd Palladium		Ag Silver		Cd Cadmium		In Indium		Sb Antimony		Sn Tin		Te Tellurium		I Iodine		Xe Xenon	
55		56		57		72		73		74		75		76		77		78		79		80		81		82		83		84		85		86	
Cs Cesium		Ba Barium		La Lanthanum		Hf Hafnium		Ta Tantalum		W Tungsten		Re Rhenium		Os Osmium		Ir Iridium		Pt Platinum		Au Gold		Hg Mercury		Tl Thallium		Pb Lead		Bi Bismuth		Po Polonium		At Astatine		Rn Radon	
87		88		89																															
Fr Francium		Ra Radium		Ac Actinium																															
						58		59		60		61		62		63		64		65		66		67		68		69		70		71			
						L Lanthanides		Ce Cerium		Pr Praseodymium		Nd Neodymium		Pm Promethium		Sm Samarium		Eu Europium		Gd Gadolinium		Tb Terbium		Dy Dysprosium		Ho Holmium		Er Erbium		Tm Thulium		Yb Ytterbium		Lu Lutetium	
						90		91		92		93		94		95		96		97		98		99		100		101		102		103			
						Ac Actinides		Th Thorium		Pa Protactinium		U Uranium		Np Neptunium		Pu Plutonium		Am Americium		Cm Curium		Bk Berkelium		Cf Californium		Es Einsteinium		Fm Fermium		Md Mendelevium		No Nobelium		Lr Lawrencium	

Analysis of human hair Results




Reference standard NIES 13 – S2 PICOFOX



Analysis of human hair Results




Reference standard NIES 13 – S4 T-STAR



Analysis of human hair Results



Hair samples results compared to literature values – S2 PICOFOX




Literature values: Tagagi et al (1986)

Analysis of human hair Results



Hair samples results compared to literature values- S2 PICOFOX




(6) Literature values: Tagagi et al (1986)

Analysis of human hair Results



Hair samples results compared to literature values- S4 T-STAR




Literature values: Tagagi et al (1986)

Analysis of human hair Results



Hair samples results compared to literature values- S4 T-STAR




(6) Literature values: Tagagi et al (1986)

Analysis of human hair Results



3-sigma LLDs



Analysis of human hair

Results & Discussion



- All elements of interest can easily be analyzed with the **S2** PICOFOX and the **S4** T·STAR, prepared as suspensions or after acid digestion
- Acid digestion gives better sensitivities but elements like Hg get lost during preparation
- The **S4** T·STAR gives distinctly better sensitivities compared to the **S2** PICOFOX
- A possible application strategy could be to use the S2 PICOFOX for a fast mobile analysis and the **S4** T·STAR for a detailed lab-based analysis
- The TXRF method generally offers a high flexibility for enhanced analytical tasks (additional nutrients, toxic metals etc.)



Part II: Application of TXRF in medical research

Selenoprotein P is the essential selenium transporter for bones†

Metallomics, 2014, 6, 1043

Nicole Pietschmann,^{ab} Eddy Rijntjes,^a Antonia Hoeg,^{ab} Mette Stoedter,^a
Ulrich Schweizer,^c Petra Seemann^d and Lutz Schomburg^{*a}

Metallomics



PAPER

[View Article Online](#)

[View Journal](#) | [View Issue](#)

Bone as matrix

Selenoprotein P is the essential selenium transporter for bones†

Metallomics, 2014,6, 1043



Nicole Pietschmann,^{ab} Eddy Rijntjes,^a Antonia Hoeg,^{ab} Mette Stoedter,^a Ulrich Schweizer,^c Petra Seemann^d and Lutz Schomburg^{*a}



Metallomics

PAPER

View Article Online
View Journal | View Issue



Se concentrations were quantified by **TXRF spectroscopy** using a **PICOFOX S2** spectrometer

(Bruker Nano GmbH, Berlin, Germany).

Briefly, bones were **decalcified** in 2.5% nitric acid and supplemented with **Gallium (Ga)** ... Alternatively, bones were treated with **proteinase K** and analysed as above.

**DOUBLE-BLIND, PLACEBO-CONTROLLED, RANDOMIZED TRIAL
OF SELENIUM IN GRAVES' HYPERTHYROIDISM**

George J Kahaly, Michaela Riedl, Jochem König, Tanja Diana, Lutz Schomburg



Se-Supplementation Trial in Humans

**DOUBLE-BLIND, PLACEBO-CONTROLLED, RANDOMIZED TRIAL
OF SELENIUM IN GRAVES' HYPERTHYROIDISM**



George J Kahaly, Michaela Riedl, Jochem König, Tanja Diana, Lutz Schomburg

	Selenium + MMI	Placebo + MMI	<i>P</i>^a
N	35	35	
Age, y, mean (SD)	44.5 (13.8)	44.5 (13.4)	1.0
Height, m, mean (SD)	1.69 (0.091)	1.69 (0.066)	0.988
Weight, kg, mean (SD)	68.8 (11.2)	71.8 (11.7)	0.270
Systolic blood pressure, mm Hg, mean (SD)	118.0 (11.5)	120.4 (14.2)	0.435

Table 2. Serological Results

	Selenium + MMI					Placebo + MMI				
	Week 0	Week 4	Week 12	Week 24	Week 36	Week 0	Week 4	Week 12	Week 24	Week 36
Selenium, µg/L	109 (25)	189 (209)	209 (261)	160 (51)	117 (43)	115 (29)	107 (34)	107 (27)	111 (27)	110 (26)
Vs baseline		<i>P</i> = 0.030	<i>P</i> = 0.035	<i>P</i> < 0.001	<i>P</i> = 0.395		<i>P</i> = 0.036	<i>P</i> = 0.024	<i>P</i> = 0.574	<i>P</i> = 0.373
Selenium vs placebo	<i>P</i> = 0.384	<i>P</i> = 0.026	<i>P</i> < 0.001	<i>P</i> < 0.001	<i>P</i> = 0.406					
SELENOP, mg/L	3.8 (0.8)	5.2 (1.1)	5.1 (1.1)	5.1 (1.3)	3.3 (0.9)	3.7 (1.1)	3.6 (1.0)	3.6 (1.0)	3.6 (1.0)	2.9 (0.8)
Vs baseline		<i>P</i> < 0.001	<i>P</i> < 0.001	<i>P</i> < 0.001	<i>P</i> = 0.008		<i>P</i> = 0.248	<i>P</i> = 0.413	<i>P</i> = 0.509	<i>P</i> < 0.001
Selenium vs placebo	<i>P</i> = 0.828	<i>P</i> < 0.001	<i>P</i> < 0.001	<i>P</i> < 0.001	<i>P</i> = 0.082					

10 µl of serum was spiked with an internal Gallium and applied to polished quartz glass carriers.
Inter-assay and intra-assay CV were below 15% as determined using a commercial human reference serum sample (Serorm; SERO AS, Billingstad, Norway).

Selenium and copper status - potential signposts for neurological remission after traumatic spinal cord injury

Julian Seelig^a, Raban Arved Heller^{a,b}, Julian Hackler^a, Patrick Haubruck^{b,c}, Arash Moghaddam^d, Bahram Biglari^e, Lutz Schomburg^{a,*}








ELSEVIER


Diagnostic Trial of serum Se and Cu in TSCI

Selenium and copper status - potential signposts for neurological remission after traumatic spinal cord injury

Julian Seelig^a, Raban Arved Heller^{a,b}, Julian Hackler^a, Patrick Haubruck^{b,c}, Arash Moghaddam^d, Bahram Biglari^e, Lutz Schomburg^{a,*}



Aminoglycoside-driven biosynthesis of selenium-deficient Selenoprotein P

Kostja Renko¹, Janine Martitz¹, Sandra Hybsier¹, Bjoern Heynisch¹, Linn Voss¹, Robert A. Everley², Steven P. Gygi², Mette Stoedter¹, Monika Wisniewska¹, Josef Köhrle¹, Vadim N. Gladyshev³ & Lutz Schomburg ¹

SCIENTIFIC REPORTS 

Analytics of recombinant or purified Protein

Aminoglycoside-driven biosynthesis of selenium-deficient Selenoprotein P

Kostja Renko¹, Janine Martitz¹, Sandra Hybsier¹, Bjoern Heynisch¹, Linn Voss¹, Robert A. Everley², Steven P. Gygi², Mette Stoedter¹, Monika Wisniewska¹, Josef Köhrle¹, Vadim N. Gladyshev³ & Lutz Schomburg¹





Figure 5

SELENOP was purified by immuno-affinity and subjected to LC-MS/MS analysis. (A) Sec in SELENOP was detected almost exclusively at the positions SEC3, SEC4 and SEC5 when cells were supplemented with selenite. (B) The pattern of amino acids inserted at the three Sec codons varied strongly when cells were grown in the presence of 100 nM selenite and 200 µg/mL G418.

Pre-diagnostic copper and zinc biomarkers and colorectal cancer risk in the European Prospective Investigation into Cancer and Nutrition cohort

Carcinogenesis, 2017, Vol. 38, No. 7, 699–707

doi:10.1093/carcin/bgx051

Advance Access publication June 1, 2017

Original Article

Nested case-control Study => predictive Biomarker


Pre-diagnostic copper and zinc biomarkers and colorectal cancer risk in the European Prospective Investigation into Cancer and Nutrition cohort

Carcinogenesis, 2017, Vol. 38, No. 7, 699–707

doi:10.1093/carcin/bgx051

Advance Access publication June 1, 2017

Original Article



Bench-top total reflection X-ray fluorescence (TXRF) spectrometer (Picofox™ S2, Bruker Nano GmbH, Berlin, Germany) was used to analyse the serum sample.



Part III: Summary & Q/A

TXRF offers an ideal analytical solution for elemental analysis in medical and biological research


- Analysis of small sample amounts in the low μl -range
- Simultaneous analysis of main- and trace elements
- Simultaneous analysis of other important samples types like buffer solutions
- Instruments can be operated in normal laboratory environments (small footprint, no external gases or cooling water necessary)
- Moderate analytical demands on laboratory staff
- Low analytical and lifetime costs


Q & A



Any Questions?

Please **type in** the questions you may have for our speakers in the **Questions Box** and click **Submit**





Innovation with Integrity

Thank you for your attention!

www.bruker.com / www.s4tstar.com