

# Relevance of Oxygen and Hydrogen Determination in Additive Manufacturing and Powder Recycling

Session 46, Powders for AM-Recycling



**Dr. Peter Papelewski \***  
**Christian Zühlke**

Bruker AXS GmbH, Elemental Analysis,  
Karlsruhe, Germany



**Max Lutter-Günther**

Fraunhofer Research Institution for  
Casting, Composite and Processing  
Technology, Augsburg, Germany

# Bruker Solutions for AM/PM Elemental and Structural Analysis



**X-Ray Fluorescence (WD-XRF)**  
Elemental composition &  
mapping



S8 TIGER

**Inert Gas Fusion –  
Mass Spectrometry (IGF-MS)**  
O, N, H & Ar detection



G6 LEONARDO  
G8 GALILEO



**Combustion Analysis**  
C and S analysis in solids




G4 ICARUS

**Electron Microscope Analyzers**  
Highest resolution Imaging  
Topographical & material  
analysis



EM ANALYZERS

**Micro Computed Tomography**  
Metrology  
Non-destructive inspection  
CAD comparison



SKYSCAN 1272

**X-Ray Diffraction (XRD)**  
Phase composition  
Texture  
Residual stress



D8 DISCOVER

**Mechanical/Tribology Test**  
Wear rates & resistance  
Friction  
Load  
Hardness



UMT TRIBOLAB

**Nano Indentation (NI)**  
Mechanical homogeneity  
Fracture  
Surface energy  
Agglomeration




TI 980 TRIBOINDENTER

**3D Optical Microscopy**  
Roughness  
Metrology  
Surface geometry  
& shape



NPFLEX

# Relevance of Oxygen & Hydrogen Introduction



- Chemical reactions on the metal particle surface play an important role during:
  - Powder Production
  - Storage
  - Conditioning for AM
  - Handling in the machine
  - Post Processing / Heat Treatment
  - Powder Recycling

⇒ **Powder Aging** (degradation, oxidation, uptake, corrosion, ...)

# Relevance of Oxygen & Hydrogen

## Effects of “gases in metals”



Laser Beam Melting is related to welding – layer by layer - and shares the same challenges

### **Oxygen**

- unwanted, parasite element
- steel: causing ageing brittleness
- Titanium: degrades mech. properties (e.g. ductility)
- bulk (inclusions) vs. surface (oxidation)
- reactive metals like Ti, Mg, Al have high oxygen affinity

### **Hydrogen**

- highly unwanted element
- causes harmful, complex embrittlement damages
- can be supplied during: welding, etching, annealing, corrosion for moisture

# Relevance of Oxygen & Hydrogen Factors for Degradation



- Powder degradation will be faster:
  - the smaller the particles are (high specific surface)
  - the higher moisture / oxygen levels are
  - the higher the temperature is
  - the longer the exposition time is




Correlation between particle size and O, H content on AlSi10Mg:

Particle Size / $\mu\text{m}$	Oxygen / ppm	Hydrogen / ppm
10 - 53	1640	68.2
20 - 63	998	53.7
40 - 100	772	38.1

# Relevance of Oxygen & Hydrogen Surface reaction with moisture





- AlSi10Mg reacts with ambient moisture
- Hydroxide layers are formed on the surface:



"Caking" on AlSi10Mg powder due to crystal connection of neighboring particles

Hydroxide layer formation on AlSi10Mg (verified by FTIR)



# Relevance of Oxygen & Hydrogen Effect on LBM Process and Quality






- New powder (A)
- Powder (B) aged: storage at high humidity and elevated temperatures
- Note: Powder (B) still shows unchanged fluidity

Sample specimen:  
A (left), B (right)

**A:** Density 99.5%

**B:** Density 96.5%





# Relevance of Oxygen & Hydrogen

## Example: Spatter Formation



### Origin


- 7.5 mg spatter material per g material fused
- Spatter end-up in powder bed



Spherical spatter, solidified in flight

### Oxygen Content [ppm]

IGF,  $n=5$



### Oxide Layer Thickness

(XPS-measurement as  $\text{Al}_2\text{O}_3$ )

Virgin powder: 27 nm

Spatter: 35 nm



# Relevance of Oxygen & Hydrogen

## How to determine O, N, H (and Ar)?






Graphite crucible on lower electrode




- Techniques like WD-XRF and OES cannot determine O, H, N or Ar
- Inert gas fusion (**IGF**) was applied in this study
  - A sample of any shape (powder, pieces, drillings, etc.) is weighed
  - The sample is fused inside a graphite crucible (up to 3000 dC) in a flow of high purity inert gas

Ti6Al4V test specimen produced by LBM



# Relevance of Oxygen & Hydrogen IGF: Temperature matters



# Relevance of Oxygen & Hydrogen

## Characteristics of IGF



### **Inert gas fusion is:**

- a volumetric method (entire specimen is analyzed)
- a relative method that requires calibration (CRM or gas dosing)  $\Rightarrow$  result traceability
- not limited in concentration ranges: high-ppb to 100%
- Fast (30s pre-cleaning of crucible, 60-90s measuring time)
- applicable to all inorganic solids (must be dry!)
- flexible by coupling a MS: Ar determination or isotopic tracers
- able to deliver kinetic information by applying heating rates (e.g. separation of oxides)

# Relevance of Oxygen & Hydrogen CRMs for AM



ISO Certified · 9001 · 17025 · 17043 · 17034

## Certificate of Analysis IARM Fe316LP-18

Additive Manufacturing Powder (-270M+16µ) Stainless Steel 316L / UNS S31603

### Certified Reference Material

Certified Values listed in wt.% with associated uncertainties

Al	0.006 ± 0.005	C	0.006 ± 0.001	Co	0.006 ± 0.003	Cr	17.9 ± 0.2
Cu	0.0031 ± 0.0005	Mn	1.56 ± 0.03	Mo	2.81 ± 0.05	N	0.081 ± 0.009
Ni	13.9 ± 0.2	O	0.043 ± 0.007	P	0.011 ± 0.003	S	0.0041 ± 0.0002
Si	0.29 ± 0.02	V	0.006 ± 0.005				



ISO Certified · 9001 · 17025 · 17043 · 17034

## Certificate of Analysis IARM Ti64P-18

Additive Manufacturing Powder (-53/+16µ) Titanium Alloy 6-4 / UNS R56400

### Certified Reference Material

Certified Values listed in wt.% with associated uncertainties

Al	6.47 ± 0.09	C	0.051 ± 0.004	Fe	0.216 ± 0.005	H	0.0018 ± 0.0007
Mn	0.011 ± 0.002	N	0.04 ± 0.02	O	0.15 ± 0.02	S	0.0014 ± 0.0006
Sn	0.008 ± 0.001	V	4.24 ± 0.05				

H	O
0.001	0.0725
0.0012	0.151
0.0013	0.154
0.0014	0.1562
0.0015	0.1593
0.0016	0.16
0.002	0.1605
0.0021	0.167
0.004	0.186
0.0018	0.15
0.0009	0.03
0.0018	0.15
0.0007	0.02
F	F



ISO Certified · 9001 · 17025 · 17043 · 17034

## Certificate of Analysis IARM Ni718P-18

Additive Manufacturing Powder (-270M+16µ) Nickel Alloy 718 / UNS N07718

### Certified Reference Material

Certified Values listed in wt.% with associated uncertainties

Al	0.49 ± 0.02	C	0.036 ± 0.003	Co	0.097 ± 0.006	Cr	19.6 ± 0.2
Cu	0.018 ± 0.004	Fe	17.0 ± 0.3	Mn	0.026 ± 0.002	Mo	3.13 ± 0.04
N	0.010 ± 0.007	Nb	4.95 ± 0.08	Ni	53.6 ± 0.3	O	0.014 ± 0.002
P	0.006 ± 0.003	S	0.0013 ± 0.0004	Si	0.036 ± 0.009	Ta	0.006 ± 0.005
Ti	1.01 ± 0.01	V	0.017 ± 0.003	W	0.010 ± 0.005		

Jernkontoret swerea | KIMAB



### Certified Reference Material

Jernkontoret

JK 47A

Iron powder

Materials supplier: Höganäs Sweden AB

CERTIFIED VALUES - Mass content in %

	O	N	C
M <sub>m</sub>	0.69	0.0062	0.370
C(95%)	0.02	0.0002	0.003

# Relevance of Oxygen & Hydrogen Ar? $\Rightarrow$ Process Gases Inclusions



- During the gas atomization, liquid metal is hit/"nebulized" by a jet of a process gas (argon or nitrogen)
- Impurities in the gas like moisture, oxygen, etc. can react
- Additionally, traces of the process gas remain in the particles:
  - Ti, Al, Mg and other reactive alloys react with nitrogen  
 $\Rightarrow$  argon is used as atomizing gas
  - **Ar** does not react but stays in the powder in form of **entrapped closed porosity**
  - Pores filled with Ar cause similar embrittlement like hydrogen but starting at much smaller levels: <100 ppb [1, 3]
  - Micro CT ( $\mu$ CT) can't tell whether pore is filled with Ar
  - Post-processing cannot remove a pore filled with gas [4]
- IGF coupled to a mass spectrometer allows ppb detection limits of argon

# Relevance of Oxygen & Hydrogen Motivation



From the AM user perspective the powder metal market is complex and non-transparent:

- Manufacturers of AM machines offer their own powder
- Premium powder manufacturers
- Smaller primary metal manufacturers starting own powder productions

Concerning the chemical composition

- Alloying elements and material grade is given and specified (and often verified multiple times – although that is invariant during the process)
- Very seldom light elements like C, S and O, H are listed
- Hydrogen and Argon are ignored (except: NADCAP)
- No quality threshold values are defined

**A lack of knowledge currently prevails regarding initial O,N,H contents of commercially available powders for powder bed fusion processes**

# Relevance of Oxygen & Hydrogen AlSi10Mg – Series I



- Five **different lots** of AlSi10Mg - from 2 years
- Same supplier/manufacturer of LBM-machine
- Identical specification (but no O,N,H given in certificate)
- Fresh powder
- $n=5$ , one STD in brackets as measure for homogeneity

Lot (year)	Oxygen / ppm	Nitrogen / ppm	Hydrogen / ppm
A (2016)	1323 (187)	55 (6)	56 (4)
B (2016)	1730 (45)	42 (14)	57 (3)
C (2017)	1074 (71)	66 (24)	56 (6)
D (2017)	1378 (100)	55 (9)	64 (3)
E (2017)	937 (18)	48 (14)	46 (3)

# Relevance of Oxygen & Hydrogen AlSi10Mg – Series II



- Identical specification - but no O,N,H given
- Five **different suppliers**, same year
- Fresh powder
- $n=5$ , one STD in brackets as measure for homogeneity

Supplier	Oxygen / ppm	Nitrogen / ppm	Hydrogen / ppm	Comment
1	1378 (100)	55 (9)	64 (3)	
2	1106 (79)	25 (9)	69 (4)	
3	2594 (119)	41 (9)	94 (6)	obtrusive in AM
4	1004 (135)	28 (9)	47 (2)	
5	1328 (43)	17 (6)	62 (5)	

Best quality AlSi10Mg  $\Rightarrow$  O:  $\leq 1000$  ppm, H:  $\sim 50$  ppm



# Relevance of Oxygen & Hydrogen

## Series III – Other Metals



- $n=5$ , one STD in brackets as measure for homogeneity

Material	Description	O / ppm	N / ppm	H / ppm	Questions...
Steel 1.7131 16MnCr5	low alloyed steel	1474 (133)	92 (3)	12 (1)	
SS316L	316L low carbon stainless steel	476 (46)	<b>112</b> (5)	10 (1)	Ar atomized?
SS316L	316L low carbon stainless steel	543 (43)	<b>675</b> (6)	11 (1)	N <sub>2</sub> atomized?
1.2709	Tool steel, highly alloyed	504 (16)	83 (2)	14 (3)	
Inconel 718		<b>358</b> (29)	<b>330</b> (6)	<b>59</b> (4)	???
Inconel 718		<b>175</b> (11)	191 (3)	<b>4.7</b> (0.2)	
Ni-Base Pwd.	comparable to Inconel	267 (19)	78 (3)	14 (1)	
18Ni300	Steel with high Ni	403 (23)	212 (1)	18 (3)	
CuCr1Zr	90%Cu, "Zr aged"	<b>2072</b> (123)	24 (3)	13 (1)	green Laser req.

# Relevance of Oxygen & Hydrogen

## Series IV: Ti-6Al-4V



### Typical values for Ti6Al4V analyzed so far (fresh powder)

- **O**: 950 – 1250 ppm (but up to 3600 ppm in specimen)  
ASTM limit grade 1: <1800 ppm
- **N**: 150 – 350 ppm
- **H**: 13 – 220 ppm
- **Ar**: 0.3 – 1.3 ppm (limit in HIP-AM: 0.1 ppm)

# Relevance of Oxygen & Hydrogen Conclusions



- Oxygen and Hydrogen are the most critical elements for quality and costs (e.g. recycling)
- To optimize the process and find the best compromise, Process- and Quality-Control of these elements is indicated
- Metal powder produced by gas atomization infiltrates process gases. Wish: The contaminants shall be specified and also tested in the final product (until experience values for AM are set)
- A high level of care concerning the used process gases and especially their impurities, like **moisture**, is required to run a robust AM process with highest quality
- Inert gas fusion is a fast, effective and readily available analytical method for process and quality control in powder metallurgical processes

## Thank you for your attention!

### References

- [1] B. Hofer, F. Bigolin, M. Hamentgen, H. Stremming, C. Zühlke, "Measuring the gas content in HIP components and impurities in the argon- and chemical reacted gas used to compacting near net shape parts and castings. Presentation of the measurement technique", Session 215, HIP 2014 conference (Stockholm, Sweden)
- [2] P. J. Mohr, B. N. Taylor, D. B. Newell; "CODATA recommended values of the fundamental physical constants - 2010"; National Institute of Standards and Technology, USA, <http://physics.nist.gov/constants>
- [3] P. Mellin, M. Östlund, W. Fredriksson, C. Pellegrini, H. Blom, O. Björnberg, A. Strondl; „Detecting Argon Trapped In Reference Samples Made By Hot Isostatic Pressing“, WorldPM 2016 (Hamburg, Germany), WP1603294934
- [4] C. Schaak, W. Tillmann, M. Schaper, M. E. Aydinöz; "Process Gas Infiltration in Inconel 718 Samples during SLM Processing", RTeJournal - Fachforum für Rapid Technologie, Vol. 2016 (2016), Germany