

Quantification of Steels and Alloys using a dual source multidetector system.

Part II: SEM-WDS adding to XRF-EDS and SEM-EDS analysis

Bruker Nano Analytics, Berlin, Germany
Webinar, April 29, 2021

Presenters



» **Andrew Menzies, PhD**
Sr. Applications Geology and Mining,
Bruker Nano Analytics, Berlin, Germany




» **Dr. Michael Abratis**
Sr. Applications Scientist WDS,
Bruker Nano Analytics, Berlin, Germany

Overview

Quantification of Steels and Alloys

Webinar structure

- **Part I** of the webinar series focused on the dual-beam sources (electron and X-ray source), and how they interact with the samples of interest to generate X-rays which are identified and quantified using EDS.
- **Part II** will now compare the earlier results with the measurements using the WDS collected on the same system.
- Multi method approach
- SEM WDS Introduction
- Samples and Methods
- Application to Steels and Alloys
- Summary and Conclusion



Webinar: Part I – available at:
<https://www.bruker.com/en/news-and-events/webinars/2021/XRF-EDS-and-SEM-EDS-Analysis.html>

Overview

Motivation for multi-method approach

- Why use different analytical techniques for steel analysis?
- Quantification of Steels and Alloys can be problematic due to the various elements of interest ranging from very low Z-numbers up to high Z-elements and the range of concentrations present from majors through to traces.
- **E-beam excitation / EDS detection:**
 - Very good low Z-element sensitivity / spatial resolution
 - Relatively high spectral background / limiting in sensitivity
- **X-ray excitation / EDS detection:**
 - Low spectral background / high sensitivity for higher Z-elements
 - Limitation in spatial resolution/ light element detection
- **E-beam excitation / WDS detection:**
 - High spectral resolution; improved P/B ratio compared to e-beam / EDS
 - Sequential analysis (one element at a time)


Excitation: X-ray and electron beam

Detectors: EDS and WDS

FEG SEM

Dual source excitation:
e-beam
X-ray beam

Dual Detectors:
EDS
WDS




What have we learned from Part I ?

Combined Analysis of Steel

Sample 32: AISI 422-205B

Element	Certified	MicroXRF	SEM-EDS	Combined
C	0.22			
N	0.05			
Al	0.01			
Si	0.37		0.34	0.33
P	0.01			
S	0.00			
Ti	0.00	0.003		0.00
V	0.26	0.279		0.26
Cr	11.72	11.084	11.37	11.32
Mn	0.68	0.797	0.87	0.75
Fe	83.70	83.243	84.55	83.20
Co	0.03	0.024	0.49	0.02
Ni	0.70	0.692	0.54	0.67
Cu	0.15	0.177		0.15
Nb	0.02	0.012		0.01
Mo	0.97	0.970	0.95	0.94



Analysis of Steels and Alloys: Combined Quantification

- If both electron and X-ray excitation are available, the benefits of both methods can be combined.


That is:

1. better light element sensitivity of electron excitation, e.g., from C to Si typically have smaller statistical error and better sensitivity,
2. better trace element sensitivity for heavy elements of X-ray excitation
 - *Thus, the results for each quantification method can be calculated separately, and then the results for the elements with better sensitivity and accuracy are used to calculate an improved combined quantification.*
 - However, it is still possible to improve the results with the addition of SEM-WDS
 - **Note:** the sample has to fulfil the requirements for both excitation types, i.e., the sample needs to be conductive and polished and has to be homogeneous


QUANTAX WDS

System Components

QUANTAX WDS: integral part of the QUANTAX family




XSense™ WD spectrometer



Signal processing unit SVE 6



ESPRIT 2



- Spectrum, P/B-acquisition in 'Spectrum' and 'Objects' mode
 - Mapping and LineScan
 - Quantification (SB, coupled quant possible)
 - Device control
- ... all integrated in the Esprit GUI

Introduction

Motivation for WDS application in steel analysis

- 01 Determine traces of low-Z elements in steel
 - Contents of Al, Si, P, S

- 02 Determine light elements in steel
 - Contents of N, C


- 03 Resolve spectral overlaps in steel
 - Peak overlaps of W-M with Si-K; Mo-L with S-K

- 04 Resolve spatial sample heterogeneities
 - Small-scale elemental variations, e.g. C

Analysis of Steels Samples

Steel & iron set I (ARMI)

- 15 Cr-Ni steels, 14 irons
- Variable major and trace element contents
- Certified compositions




Steel set II (ACX)

- 15 Cr-Ni steels
- Variable major and trace element contents
- Certified compositions



Analysis of Steels Samples



➤ Major Elements:

Cr, Fe, Mn, Ni

➤ Minor and Trace Elements:


C, N, Al, Si, P, S, Ti, V, Co, Cu, Nb, Mo, Sn, W

Element	Minimum(%)	Maximum(%)	Range(%)
C	0.02	1.02	1.00
N	0.01	0.33	0.32
Al	0.00	1.16	1.16
Si	0.27	1.38	1.11
P	0.01	0.04	0.03
S	0.00	0.29	0.29
Ti	0.00	0.63	0.63
V	0.02	0.26	0.24
Cr	11.72	23.60	11.88
Mn	0.35	9.31	8.96
Fe	41.29	86.23	44.95
Co	0.02	0.18	0.17
Ni	0.11	35.84	35.73
Cu	0.03	0.47	0.44
Nb	0.00	0.72	0.72
Mo	0.06	2.30	2.25
Sn	0.00	0.01	0.01
W	0.01	1.10	1.09
Total	100.00	100.00	0.00

WDS methodology

Major elements in steel


- WDS designed for highest efficiency in low X-ray energy region
- Fe, Cr, Ni determination on L-lines
- Energy selective technique → sequential acquisitions
- Acquisitions on reference material necessary
- Major elements equally covered by EDS, micro-XRF
- Possible but not required by WDS



WDS methodology

Trace elements in steel


- Better peak resolution
- High peak to background ratios
- Low detection limits
- Especially for lines $< \sim 2$ keV (complements μ XRF)
- Time effort rewarded by high precision and accuracy



WDS methodology

Two different acquisition modes


Energy range scans



Suitable for qualitative determination

= intuitive, convincing presentation

Peak to background measurements



Suitable for quantification


= faster acquisition, better statistics

Methodology

Standard-based quantification

- Determine net peak intensities of „unknown“ sample
- Compare to reference (standard) net peak intensities
- Acquire at identical conditions (geometry, acceleration voltage)
- Record beam current value
- Correct for matrix effects (ZAF)

$$\frac{C_i^S}{C_i^{Std}} \propto \frac{N_i^S}{N_i^{Std}} \frac{(k_Z k_A k_F)_{i,S}}{(k_Z k_A k_F)_{i,Std}}$$






Trace element determination by WDS

Silicon (Si-K α)

Analytical
Conditions Point
Analysis:

20 kV, 70 nA,,
under vacuum,
WD 15 mm,
300 seconds



Si Range:
0.19 – 1.99 wt.%
83% < 0.6 wt.%

WDS Results:



Minor to trace
element
concentrations

Accurate and precise

Trace element determination by WDS

Silicon (Si-K α)

Focus on lowest concentrations (0.2 – 0.7 wt.%)



Peak resolution by WDS


Tungsten (W-M α) and Silicon (Si-K α)

WDS Results:

Determination of Si and W

Comparing WDS and EDS

Peak resolution vs. peak deconvolution



Si: 0.37 wt.%

W: 1.1 wt.%

Δ Si-K α – W-M α : 35 eV




Trace element determination by WDS

Phosphorous (P-K α)

WDS Results:

Trace element concentrations

Still well matching accepted concentration values



Range:

P: 140 – 400 ppm




Trace element determination by WDS

Molybdenum (Mo-L α)

WDS Results:

Precise and accurate

Wide range of concentrations



Range:

Mo: 0.02 – 3.11 wt.%

77% < 0.52 wt.%



Peak resolution by WDS


Molybdenum (Mo-L α) and Sulfur (S-K α)

WDS Results:

Determination of Mo and S

Comparing WDS and EDS

Peak resolution vs. peak deconvolution



Mo: 0.08 wt.%

S: 0.29 wt.%


Δ S-K α - Mo-L α : 14 eV



Trace element determination by WDS

Sulfur (S-K α)

WDS Results:
Most match
accepted
values
Some deviate
due to high Mo



Range:
S: 4 – 2900 ppm
86% < 300 ppm




Light element determination by WDS

Nitrogen (N-K α) with BRML60

Analytical
Conditions
Point Analysis:

4 kV, 20 nA,,
under vacuum,
WD 15 mm,
60 seconds



Range:

N: 0.01 – 0.33 wt.%
80% < 600 ppm

BRML60
low sensitivity

25 cps/nA
P/B: 36
@N-K α , 10kV, BN




Light element determination by WDS

Nitrogen (N-K α) with BRML80

WDS Results:

BRML80 option
high sensitivity
precise
but offset
background



Range:

N: 0.01 – 0.33 wt.%

80% < 600 ppm

BRML80

250 cps/nA

P/B: 39

@N-K α , 10kV, BN

Background:

Sc-L α @ 395 eV




Light element determination by WDS

Nitrogen (N-K α) with BRML80

WDS Results:

BRML80
corrected
accurate and
precise



Range:

N: 0.01 – 0.33 wt.%

80% < 600 ppm

BRML80
250 cps/nA
P/B: 39
@N-K α , 10kV, BN

Background corrected




Light element determination by WDS

Nitrogen (N-K α) with BRML80

WDS Results:

High sensitivity
for light
elements

Precise and
accurate down
to 100 ppm



Range:

N: 0.01 – 0.33 wt.%

80% < 600 ppm

BRML80
Offset corrected




Light element determination by WDS

Carbon (C-K α)

WDS Results:

High C samples match lower limit currently at 0.3 wt.%



Range:

C: 0.002 – 1.04 wt.%
69% < 0.3 wt.%


Limitation for carbon determination

Not due to WDS sensitivity

But due to system blank

Light element determination by WDS


Carbon contamination during sample analysis




Carbon
contamination
on a steel
sample surface

Light element determination by WDS

Carbon contamination over time




Increasing carbon deposition on sample surface during measurement




Light element determination by WDS

Heterogeneous carbon distribution

Carbide bearing steel



Dual phase steel





Methodology

Combined analysis process

Element	SEM-EDS
	initial
C	
N	
Al	
Si	
P	
S	
Ti	
V	
Cr	
Mn	
Fe	
Co	
Ni	
Cu	
Nb	
Mo	



Methodology

Combined analysis process


Element	SEM-EDS	MicroXRF	SEM-EDS
	initial	iterative	refined
C		Blue	
N		Blue	
Al	Diagonal lines	Blue	Diagonal lines
Si	Light Green	Blue	Light Green
P		Blue	
S		Blue	
Ti	Diagonal lines	Light Green	
V	Diagonal lines	Light Green	
Cr	Light Green	Light Green	Light Green
Mn	Light Green	Light Green	Light Green
Fe	Light Green	Light Green	Light Green
Co	Diagonal lines	Light Green	Light Blue
Ni	Light Green	Light Green	Light Green
Cu	Diagonal lines	Light Green	Light Blue
Nb	Diagonal lines	Light Green	Light Blue
Mo	Diagonal lines	Light Green	






Methodology

Combined analysis process



Methodology


Combined analysis process





Methodology

Combined analysis process





Methodologies

Combined analysis process


Element	SEM-EDS	MicroXRF	SEM-EDS	SEM-WDS	SEM-EDS	SEM-WDS	MicroXRF	SEM-EDS	Element	Final
	initial	iterative	refined	combined	combined	refined	refined	refined		combined
C				≤ 0.36		≤ 0.36			C	≤ 0.36
N				0.02		0.02			N	0.02
Al				0.00		0.00			Al	0.00
Si	0.36		0.36	0.328	0.36	0.328		0.36	Si	0.35
P				0.029		0.029			P	0.03
S				0.013		0.013			S	0.01
Ti		0.002					0.002		Ti	0.002
V		0.059					0.059		V	0.058
Cr	17.09	17.496	17.09		17.09		17.496	17.09	Cr	17.41
Mn	1.92	1.788	1.92		1.92		1.788	1.92	Mn	1.74
Fe	69.33	70.713	69.33		69.33		70.713	69.33	Fe	70.11
Co	0.60	0.154	0.60		0.60		0.154	0.60	Co	0.14
Ni	8.86	9.012	8.86		8.86		9.012	8.86	Ni	9.00
Cu	0.50	0.501	0.50	0.458	0.50	0.458	0.501	0.50	Cu	0.49
Nb	0.74	0.695	0.74		0.74		0.695	0.74	Nb	0.67
Mo	0.44	0.448	0.44	0.393	0.44	0.393	0.448	0.44	Mo	0.42

Analysis of Steels and Alloys

Combined Analysis


Sample 32: AISI 422-205B

Element	Certified	SEM-WDS	MicroXRF	SEM-EDS	Combined
C	0.22	≤ 0.34			≤ 0.34
N	0.05	0.05			0.05
Al	0.01	0.004			0.00
Si	0.37	0.339		0.34	0.34
P	0.01	0.018			0.02
S	0.00	0.029			0.00
Ti	0.00		0.003		0.00
V	0.26		0.279		0.26
Cr	11.72		11.084	11.37	11.32
Mn	0.68		0.797	0.87	0.75
Fe	83.70		83.243	84.55	83.20
Co	0.03		0.024	0.49	0.02
Ni	0.70		0.692	0.54	0.67
Cu	0.15	0.162	0.177		0.17
Nb	0.02		0.012		0.01
Mo	0.97	0.919	0.970	0.95	0.97



SEM-WDS

Low-Z elements



SEM-XRF

High-Z elements

SEM-EDS detection



Analysis of Steels and Alloys

Combined Analysis

Sample 28: AISI 347-8D

Element	Certified	SEM-WDS	MicroXRF	SEM-EDS	Combined
C	0.05	≤ 0.36			≤ 0.36
N	0.02	0.02			0.02
Al	0.00	0.00			0.00
Si	0.36	0.328		0.36	0.35
P	0.03	0.029			0.03
S	0.03	0.013			0.01
Ti	0.00		0.002		0.002
V	0.06		0.059		0.058
Cr	17.30		17.496	17.09	17.41
Mn	1.76		1.788	1.92	1.74
Fe	69.33		70.713	69.33	70.11
Co	0.14		0.154	0.60	0.14
Ni	9.19		9.012	8.86	9.00
Cu	0.47	0.458	0.501	0.50	0.49
Nb	0.72		0.695	0.74	0.67
Mo	0.44	0.393	0.448	0.44	0.42

SEM-WDS
Low-Z elements

SEM-XRF
High-Z elements

SEM-EDS
detection

Summary and Conclusions

QUANTAX WDS Benefits

High spectral resolution

- Pathological EDS peak overlaps can be resolved (e.g. Si-W, Mo-S)
- High peak to background ratio leads to low detection limits (e.g., Al, Si, P)

High performance for low X-ray energies

- Analysis of light elements (e.g. C, N)
- Trace contents for low-Z elements (e.g., Al, Si, P)

High spatial resolution

- Allows analysis of sub- μm structures (e.g., carbides)
- Full performance also at low kV (e.g., 2-10 kV)

Summary and Conclusions: Analysis of Steels and Alloys

X-ray and Electron Excitation can work in combination to provide improved quantitative results, using the benefits of each. Specifically:

- Electron excitation is preferable for light elements, e.g. C to Si
- X-ray excitation is preferable for heavy elements and trace concentrations

Samples can be analysed and quantified either as spot (point) analyses or from hypermaps.

Benefits of each analytical method can be utilised. For example:

- **Micro-XRF:** Sample Preparation is minimal for micro-XRF
 - No carbon-coating, No polishing
- **Electron Beam:** High resolution for detecting small inclusions
- **WDS:** High spectral and spatial resolution, high performance for low X-ray energies



Thank you!

Michael Abratis
michael.abratis@bruker.com

More Information

For more information, please contact us:

Bruker Nano GmbH

info.bna@bruker.com

and

If you want to learn more about practical micro-XRF or SEM-EDS / WDS analysis including sample, measurement setup, and evaluation, our latest video series is available via the Bruker website and youtube

Product Videos



Part I - Introduction to micro-XRF and the Rapid Stage on a SEM



Part II - Loading a sample and performing a measurement



Part III - Analyzing a measured dataset

Rapid Stage




The Rapid Stage can be mounted on top of the SEM stage for fast mapping over large sample areas.



Questions and Answers

Are There Any Questions?

Please type in the questions you might have
in the Q&A box and press *Send*.



Innovation with Integrity