P162: PASEF-PRM LIVE on the timsTOF Pro

He Zhu¹, Scott B. Ficarro¹, William M. Alexander¹, Laura E. Fleming¹, Wai Cheung Chan¹, Shourjo Ghose², Matthew Willetts², Jens Decker³, Sven Brehmer³, Gary Kruppa⁴, Jarrod A. Marto^{1,*}

¹ Dana-Farber Cancer Institute, Brigham and Women's Hospital and Harvard Medical School, Boston, USA ² Bruker Daltonics Inc, Billerica, USA ³ Bruker Daltonik GmbH, Bremen, Germany ⁴ Bruker S.R.O., District Brno-City, Czech Republic

Dana-Farber HUPO Connect 2020 Virtual Congress **Cancer Institute** * email: Jarrod Marto@dfci.harvard.edu

BRIGHAM AND



BLAIJ WOMEN'S HOSPITAL



advanced center for Parkinson's disease research

Overview


Overall objective

 Develop PASEF-PRM LIVE, a framework that implements on-the-fly adjustment of the detection window to maximize target peptide coverage while responding to changes in target elution times.

General development workflow

- Initial development based on iRT spiked CSF biofluids
- Evaluate impact of real time detection window adjustment algorithms
- Evaluate reproducibility
- Explore current limits for total # peptides in 1-hr gradient
- □ Application of PRM-LIVE in drug discovery
 - Activity-based protein profiling (ABPP) for rapid selectivity analysis of inhibitors targeting deubiquitinating enzymes

Real-time RT correction algorithm



Algorithm:

- Time warping: **TW** (when ≤ 2 std peps detected)
- Time warping linear regression:
 TWLR (when ≥ 3 std peps detected)

Creating target files:


- PASEF run for CSF peptides spiked with iRTs
- Target peptides: m/z; charge; RT; 1/k0 range
- iRT peptides: m/z; RT; top6 frag ions with relative intensity

PRM LIVE:

- 112 target peptides @ 45min PASEF
- 30, 45, 60min gradients

Number of quantified peptides

gradient	#peps quantified
30min	104
45min	100
60min	103




Real-time RT correction algorithm

Quant Correlation

RT Correlation

Gradients: 30min, 45min and 60 min

Accuracy of algorithm TW/TWLR: RT-PRM-LIVE vs. RT-PASEF




Good quant correlation of targeted peptides between 30min, 45min and 60min gradients. RT of targets in PRM-LIVE has a good alignment with RT in the corresponding PASEF runs, indicating that the current real-time RT correction algorithm can accurately predict RT of targeted peptides.

Note: X annotates outliers suffering from peptides stability issue and variability of targeting window (under investigation).

Reproducibility

Quant correlation: replicates with different gradients

- 108 target CSF peptides @ 45 min PASEF
- 30, 45, 60min grads with 2 replicates (every 5.5h/replicate)



Current maximum targeted peptide# from CSF biofluids

RT tolerances


Quantified peptides # under different

20

24

Targets: 1731 CSF peps @ 60 min PASEF


RT variance \leq 1min from 3 replicates



With 30s RT tolerance, around 1400 peptides were quantified.

1352 quantified peptides were overlapped between two replicates, among which 1054 peptides (cv ≤ 0.05) showed good quant correlation.


Targeted DUBome analysis by high throughput Zip-Chip-timsTOF




Targeted protein degradation by inhibiting deubiquitination (DUB)

Human DUBome: ~ 100 enzymes

Specific DUB inhibitor screening by multiplexed ABPP



High throughput analysis






57 DUBs @ 2h on LC-timsTOF

? #DUBs @ 15 min on Zip-Chip-timsTOF

Targeted DUBome analysis by high throughput Zip-Chip-timsTOF



Summary

- Real time RT adjusting algorithm of PRM-LIVE can predict the correct RT of targeted peptides across different gradients.
- Preliminary results shows that PRM-LIVE has a good quant reproducibility of peptide targets, and PRM-LIVE is promising to target around 1,000 peptides from CSF biofluids at a 60 min LC-gradient.
- A promising initial application to activity-based protein profiling for rapid selectivity analysis of small molecular targeting deubiquitinating enzymes (DUBs); PRM-LIVE integrated with Zip-Chip separation can efficiently target peptides from 33 DUBs in a 3 min elution window (a total of 15 min run).

Acknowledgment



<u>DFCI</u>

Jarrod Marto Scott B. Ficarro Wai Cheung Chan William M. Alexander Laura E. Fleming Sharifzadeh Shabnam Joseph Card Isidoro Tavares

BRIGHAM AND

WOMEN'S HOSPITAL

Bruker Gary Kruppa Jens Decker Sven Brehmer Michael Krause Stefan Harsdorf Shourjo Ghose Matt Willetts Paul Speir

advanced center for

Parkinson's disease

research

Funding from the Michael J. Fox Foundation and the Tyler Family Foundation

BLA



