

Product Note R51-02/22

1850507 BRAVO Tripod Mount Adaptor



The BRAVO tripod mount adaptor is designed to mount the handheld Raman spectrometer BRAVO on a standard camera tripod.

The adaptor including its 120 mm z-shift allows to align the BRAVO precisely in front of objects such as art works to offer most stable measurement conditions.



Specifications

1850507 BRAVO Tripod Mount Adaptor

Allows forward/backward movement of max. 120 mm of the BRAVO for fine adjustment of the working distance.

1015487 Tripod

- Material: Aluminum and steel
- Maximum load: 12 kg
- Weight: 6.15 kg
- Min. height: 44 cm
- Max. height: 217 cm
- Closed length: 94 cm
- Inclusive: Carrying strap

1015488 Geared Head (optional)

Allows precise forward and lateral tilting and rotation of the spectrometer mounted on the tripod.

- Material: Aluminum
- Max front tilt: $\pm 30^\circ$
- Max side tilt: $\pm 7.5^\circ$
- Panoramic rotation: 360°
- Maximum load: 10 kg
- Weight: 2.6 kg

Article	Price List No.
BRAVO Tripod Mount Adaptor	1850507
Tripod	1015487
Geared Head (optional)	1015488