

RUO



## Fungiplex<sup>®</sup> Candida Auris

● Real-Time PCR

# Fungiplex Candida Auris

Fungiplex Candida Auris is a real-time PCR test, which specifically targets *Candida auris*; it provides rapid results with high sensitivity and specificity without isolation of the organism.

## Specific detection of *Candida auris*

*C. auris* is an emerging global health threat.<sup>1</sup> Since its isolation in 2009, there have been cases reported in at least 15 countries over 5 continents. Transmissible in healthcare settings, *C. auris* can cause outbreaks and is often multi-drug resistant. While some laboratory methods misidentify the organism, Bruker's MALDI Biotyper® is authorized by the FDA for use in the unambiguous identification of *Candida auris*, upon isolation of the organism. PCR also shows promise as a sensitive and specific tool for *C. auris* detection. By providing rapid results without isolation of the organism, *Candida auris* PCR could prove advantageous for infection control and patient management.

The Fungiplex Candida Auris assay is designed to rapidly detect the presence of *C. auris* DNA extracted directly from samples, without culture. Fungiplex Candida Auris is specific for *C. auris* and does not cross-react with other, closely-related, *Candida* species.

<sup>1</sup> A. Chowdhary et al. PLoS Pathog. (2017) 13(5): e1006290 DOI: <https://doi.org/10.1371/journal.ppat.1006290>

Target <sup>a</sup>	Limit of detection <sup>b</sup>	Thermocyclers (validated)	Time to result
<i>Candida auris</i>	10 input copies (gDNA)	√ ABI 7500 √ ABI Quantstudio 5 √ Bio-Rad CFX √ Hain FluoroCycler XT √ Mic qPCR Cyclers √ Qiagen RotorGene √ Roche LightCycler 480 II	< 2 h from extraction (1 - 50 samples)

<sup>a</sup> 28 different *C. auris* strains detected in triplicate. Does not cross-react with *C. haemulonii*/*C. haemulonis*, *C. duobushaemulonii*, *C. dubliniensis*, *C. lipolytica*, *C. lusitanae*, *C. albicans*, *C. glabrata*, *C. krusei*, *C. sake*, *C. parapsilosis*

<sup>b</sup> LOD refers to 100 % detection

## Fungiplex Candida Auris

- has been designed in a **user-friendly** format that is compatible with **existing laboratory equipment** and produces results from extracted nucleic acid in **less than 2 hours**.
- offers the potential to rapidly detect *C. auris* without isolation of the organism.

**Order Information: Part-No. 1860224** Fungiplex Candida Auris RUO PCR Kit

Real-time PCR Kit for the detection of *Candida auris* from extracted DNA.  
The kit provides all necessary reagents for 50 tests.

Fungiplex is a registered trademark of the Bruker group of companies.

For Research Use Only. Not for use in clinical diagnostic procedures.  
Please contact your local representative for availability in your country.



**As of May 2021, Bruker Daltonik GmbH is now Bruker Daltonics GmbH & Co. KG.**

●  **Bruker Daltonics GmbH & Co. KG** **Bruker Scientific LLC**

Bremen · Germany  
Phone +49 (0) 421-2205-0

Billerica, MA · USA  
Phone +1 (978) 663-3660

[info.md@bruker.com](mailto:info.md@bruker.com) - [www.bruker.com/microbiology](http://www.bruker.com/microbiology)