



Bruker-Hain Diagnostics

# FluoroType<sup>®</sup> MTBDR V2.0

## NEW – WHO recommendation!

Multidrug-resistant tuberculosis (MDR-TB) is an alarming and ongoing global issue which requires fast intervention in order to prevent further spread of the disease. We are committed to fighting TB worldwide and are excited that the WHO approves and recommends our **FluoroType<sup>®</sup> MTBDR V2.0 particularly for its comprehensive detection of multi-drug resistances:**

### FluoroType<sup>®</sup> MTBDR VER2.0

Procurement is now eligible using The Global Fund as it has passed their Quality Assurance Policy for Diagnostics Products

- Recognised by WHO Guideline Development Group to be highly accurate when using respiratory samples from people with signs and symptoms of pulmonary TB
- Based on the innovative LiquidArray<sup>®</sup> technology enabling highly complex results to evaluate more than 50 drug resistance gene mutations

Please contact your local representative for further details and availability in your country.

### Order Information:

FluoroType<sup>®</sup> MTBDR VER2.0  
FluoroCycler<sup>®</sup> XT

**Part No. H-62896**  
**Part No. H-2017**

Please contact your local representative for availability in your country. Not for sale in the USA.



Hardwiesenstrasse 1  
72147 Nehren - Germany  
Phone +49 (0) 7473-9451-0



[info.mdx.de@bruker.com](mailto:info.mdx.de@bruker.com) - [www.hain-lifescience.de](http://www.hain-lifescience.de)